

Phospho-IGF-IR (Y1165/Y1166) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10517	SYNONYMS
SIZE	100 ug	IGF1R; Insulin-like growth factor 1 receptor; Insulin-like
HOST	Rabbit	growth factor I receptor; IGF-I receptor; CD antigen CD221

PRODUCT INFORMATION

GENE NAME	IGF1R
PROTEIN NAME	Insulin-like growth factor 1 receptor
IMMUNOGEN	Synthesized peptide derived from human IGF-IR around the phosphorylation site of Y1165/Y1166.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.