

Phospho-I κ B- α (S32/S36) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10527	SYNONYMS
SIZE	100 ug	NFKBIA; IKBA; MAD3; NFKBI; NF-kappa-B inhibitor alpha; I-kappa-B-alpha; IkB-alpha; IkappaBalpha; Major histocompatibility complex enhancer-binding protein
HOST	Rabbit	MAD3

PRODUCT INFORMATION

GENE NAME	NFKBIA
PROTEIN NAME	NF-kappa-B inhibitor alpha
IMMUNOGEN	Synthesized peptide derived from human I κ B- α around the phosphorylation site of S32/S36.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.