

Phospho-Met (Y1349) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10543	SYNONYMS
SIZE	100 ug	MET; Hepatocyte growth factor receptor; HGF receptor;
HOST	Rabbit	HGF/SF receptor; Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met

PRODUCT INFORMATION

GENE NAME	MET
PROTEIN NAME	Hepatocyte growth factor receptor
IMMUNOGEN	Synthesized peptide derived from human Met around the phosphorylation site of Y1349.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.