

Phospho-Op18 (S38) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10558	SYNONYMS
SIZE	100 ug	STMN1; C1orf215; LAP18; OP18; Stathmin; Leukemia-associated phosphoprotein p18; Metablastin; Oncoprotein
HOST	Rabbit	18; Op18; Phosphoprotein p19; pp19; Prosolin; Protein Pr22; pp17

PRODUCT INFORMATION

GENE NAME	STMN1
PROTEIN NAME	Stathmin
IMMUNOGEN	Synthesized peptide derived from human Op18 around the phosphorylation site of S38.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.