

Phospho-p53 (S20) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10564	SYNONYMS
SIZE	100 ug	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-
HOST	Rabbit	13; Phosphoprotein p53; Tumor suppressor p53

PRODUCT INFORMATION

GENE NAME	TP53
PROTEIN NAME	Cellular tumor antigen p53
IMMUNOGEN	Synthesized peptide derived from human p53 around the phosphorylation site of S20.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.