

Phospho-PLC γ 1 (Y783) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10578	SYNONYMS
SIZE	100 ug	PLCG1; PLC1; 1-phosphatidylinositol 4; 5-bisphosphate phosphodiesterase gamma-1; PLC-148; Phosphoinositide phospholipase C-gamma-1; Phospholipase C-II; PLC-II; Phospholipase C-gamma-1; PLC-gamma-1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PLCG1
PROTEIN NAME	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1
IMMUNOGEN	Synthesized peptide derived from human PLC γ 1 around the phosphorylation site of Y783.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.