

Phospho-Sp1 (T739) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10585	SYNONYMS
SIZE	100 ug	SP1; TSFP1; Transcription factor Sp1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	SP1
PROTEIN NAME	Transcription factor Sp1
IMMUNOGEN	Synthesized peptide derived from human Sp1 around the phosphorylation site of T739.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.