

Phospho-VASP (S157) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10596	SYNONYMS
SIZE	100 ug	VASP; Vasodilator-stimulated phosphoprotein; VASP
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	VASP
PROTEIN NAME	Vasodilator-stimulated phosphoprotein
IMMUNOGEN	Synthesized peptide derived from human VASP around the phosphorylation site of S157.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.