

Phospho-Ephrin-B1/2 (Y330) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10601	SYNONYMS
SIZE	100 ug	EFNB1; EFL3; EPLG2; LERK2; Ephrin-B1; EFL-3;
HOST	Rabbit	ELK ligand; ELK-L; EPH-related receptor tyrosine kinase ligand 2; LERK-2; EFNB2; EPLG5; HTKL; LERK5; Ephrin-B2; EPH-related receptor tyrosine kinase ligand 5; LERK-5; HTK ligand; HTK-L

PRODUCT INFORMATION

GENE NAME	EFNB1/EFNB2
PROTEIN NAME	Ephrin-B1/Ephrin-B2
IMMUNOGEN	Synthesized peptide derived from human Ephrin-B1/2 around the phosphorylation site of Y330.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.