

Phospho-Crk II (Y221) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10630	SYNONYMS
SIZE	100 ug	CRK; Adapter molecule crk; Proto-oncogene c-Crk; p38
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CRK
PROTEIN NAME	Adapter molecule crk
IMMUNOGEN	Synthesized peptide derived from human Crk II around the phosphorylation site of Y221.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.