

Phospho-GABAA R β 1 (S434) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10634	SYNONYMS
SIZE	100 ug	GABRB1; Gamma-aminobutyric acid receptor subunit beta-1; GABA(A) receptor subunit beta-1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	GABRB1
PROTEIN NAME	Gamma-aminobutyric acid receptor subunit beta-1
IMMUNOGEN	Synthesized peptide derived from human GABAA angle: R β 1 around the phosphorylation site of S434.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.