

## Phospho-Shc (Y349) Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10635	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	SHC1; SHC; SHCA; SHC-transforming protein 1; SHC-transforming protein 3; SHC-transforming protein A; Src
<b>HOST</b>	Rabbit	homology 2 domain-containing-transforming protein C1; SH2 domain protein C1

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	SHC1
<b>PROTEIN NAME</b>	SHC-transforming protein 1
<b>IMMUNOGEN</b>	Synthesized peptide derived from human Shc around the phosphorylation site of Y349.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB,IHC,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**