

Phospho-Cytokeratin 18 (S33) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10690	SYNONYMS
SIZE	100 ug	KRT18; CYK18; PIG46; Keratin; type I cytoskeletal 18;
HOST	Rabbit	Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18

PRODUCT INFORMATION

GENE NAME	KRT18
PROTEIN NAME	Keratin type I cytoskeletal 18
IMMUNOGEN	Synthesized peptide derived from human Cytokeratin 18 around the phosphorylation site of S33.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.