

3pK Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10729	SYNONYMS
SIZE	100 ug	MAPKAPK3; MAP kinase-activated protein kinase 3;
HOST	Rabbit	MAPK-activated protein kinase 3; MAPKAP kinase 3; MAPKAP-K3; MAPKAPK-3; MK-3; Chromosome 3p kinase; 3pK

PRODUCT INFORMATION

GENE NAME	MAPKAPK3
PROTEIN NAME	MAP kinase-activated protein kinase 3
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human 3pK. at AA range: 270-350

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.