

4.1G Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10731	SYNONYMS
SIZE	100 ug	EPB41L2; Band 4.1-like protein 2; Generally expressed
HOST	Rabbit	protein 4.1; 4.1G

PRODUCT INFORMATION

GENE NAME	EPB41L2
PROTEIN NAME	Band 4.1-like protein 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human 4.1G. at AA range: 530-610

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.