

## Abl1/2 Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10758	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	ABL1; ABL; JTK7; Tyrosine-protein kinase ABL1;
<b>HOST</b>	Rabbit	Abelson murine leukemia viral oncogene homolog 1; Abelson tyrosine-protein kinase 1; Proto-oncogene c- Abl; p150; ABL2; ABLL; ARG; Abelson tyrosine- protein kinase 2; Abelson murine leukemia vira

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	ABL1/ABL2
<b>PROTEIN NAME</b>	Tyrosine-protein kinase ABL1/2
<b>IMMUNOGEN</b>	Synthesized peptide derived from human Abl1/2 around the non-phosphorylation site of Y393/439.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse
<b>TESTED APPLICATIONS</b>	WB,IHC,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**