

ADAR2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10782	SYNONYMS
SIZE	100 ug	ADARB1; ADAR2; DRADA2; RED1; Double-stranded
HOST	Rabbit	RNA-specific editase 1; RNA-editing deaminase 1; RNA-editing enzyme 1; dsRNA adenosine deaminase

PRODUCT INFORMATION

GENE NAME	ADARB1
PROTEIN NAME	Double-stranded RNA-specific editase 1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human ADAR2. at AA range: 450-530

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.