

AIFL Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10796	SYNONYMS
SIZE	100 ug	AIFM3; AIFL; Apoptosis-inducing factor 3; Apoptosis-inducing factor-like protein
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	AIFM3
PROTEIN NAME	Apoptosis-inducing factor 3
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human AIFL. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.