

CLK4 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11185	SYNONYMS
SIZE	100 ug	CLK4; Dual specificity protein kinase CLK4; CDC-like
HOST	Rabbit	kinase 4

PRODUCT INFORMATION

GENE NAME	CLK4
PROTEIN NAME	Dual specificity protein kinase CLK4
IMMUNOGEN	Synthesized peptide derived from the Internal region of human CLK4. at AA range: 70-150

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.