

CLP36 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11187	SYNONYMS
SIZE	100 ug	PDLIM1; CLIM1; CLP36; PDZ and LIM domain protein
HOST	Rabbit	1; C-terminal LIM domain protein 1; Elfin; LIM domain protein CLP-36

PRODUCT INFORMATION

GENE NAME	PDLIM1
PROTEIN NAME	PDZ and LIM domain protein 1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human CLP36. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.