

CMTM3 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11190	SYNONYMS
SIZE	100 ug	CMTM3; CKLFSF3; CKLF-like MARVEL
HOST	Rabbit	transmembrane domain-containing protein 3; Chemokine-like factor superfamily member 3

PRODUCT INFORMATION

GENE NAME	CMTM3
PROTEIN NAME	CKLF-like MARVEL transmembrane domain-containing protein 3
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human CMTM3. at AA range: 1-80

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.