

CNT2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11198	SYNONYMS
SIZE	100 ug	SLC28A2; CNT2; Sodium/nucleoside cotransporter 2;
HOST	Rabbit	Concentrative nucleoside transporter 2; CNT 2; hCNT2; Na(+)/nucleoside cotransporter 2; Sodium-coupled nucleoside transporter 2; Sodium/purine nucleoside co- transporter; SPNT; Solute carrier

PRODUCT INFORMATION

GENE NAME	SLC28A2
PROTEIN NAME	Sodium/nucleoside cotransporter 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human CNT2. at AA range: 340-420

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.