

Cofilin Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11200	SYNONYMS
SIZE	100 ug	CFL1; CFL; Cofilin-1; 18 kDa phosphoprotein; p18;
HOST	Rabbit	Cofilin; non-muscle isoform; CFL2; Cofilin-2; Cofilin, muscle isoform

PRODUCT INFORMATION

GENE NAME	CFL1/CFL2
PROTEIN NAME	Cofilin-1/2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Cofilin. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.