

COL11A2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11201	SYNONYMS
SIZE	100 ug	COL11A2; Collagen alpha-2(XI) chain
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	COL11A2
PROTEIN NAME	Collagen alpha-2(XI) chain
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human COL11A2. at AA range: 1180-1260

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.