

COL1A2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11202	SYNONYMS
SIZE	100 ug	COL1A2; Collagen alpha-2(I) chain; Alpha-2 type I collagen
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	COL1A2
PROTEIN NAME	Collagen alpha-2(I) chain
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human COL1A2. at AA range: 1-80

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.