

COL4A6 Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|------------------------------------|
| CATALOG NO. | PTX11206 | SYNONYMS |
| SIZE | 100 ug | COL4A6; Collagen alpha-6(IV) chain |
| HOST | Rabbit | |

PRODUCT INFORMATION

| | |
|---------------------|---|
| GENE NAME | COL4A6 |
| PROTEIN NAME | Collagen alpha-6(IV) chain |
| IMMUNOGEN | Synthesized peptide derived from the C-terminal region of human COL4A6. at AA range: 1170-1250 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human |
| TESTED APPLICATIONS | WB,IHC,IF,ELISA |
| RECOMMENDED DILUTION | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.