

Connexin 25 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11209	SYNONYMS
SIZE	100 ug	GJB7; CX25; Gap junction beta-7 protein; Connexin-25;
HOST	Rabbit	Cx25

PRODUCT INFORMATION

GENE NAME	GJB7
PROTEIN NAME	Gap junction beta-7 protein
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Connexin 25. at AA range: 40-120

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.