

## CPN cat Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX11223	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	CPN1; ACBP; Carboxypeptidase N catalytic chain; CPN;
<b>HOST</b>	Rabbit	Anaphylatoxin inactivator; Arginine carboxypeptidase; Carboxypeptidase N polypeptide 1; Carboxypeptidase N small subunit; Kininase-1; Lysine carboxypeptidase; Plasma carboxypeptidase B; S

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	CPN1
<b>PROTEIN NAME</b>	Carboxypeptidase N catalytic chain
<b>IMMUNOGEN</b>	Synthesized peptide derived from the C-terminal region of human CPN cat. at AA range: 380-460

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human
<b>TESTED APPLICATIONS</b>	WB,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**