

CREB-2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11229	SYNONYMS
SIZE	100 ug	ATF4; CREB2; TXREB; Cyclic AMP-dependent transcription factor ATF-4; cAMP-dependent transcription factor ATF-4; Activating transcription factor 4; Cyclic AMP-responsive element-binding protein 2; CREB-2; cAMP-responsive element-binding prot
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	ATF4
PROTEIN NAME	Cyclic AMP-dependent transcription factor ATF-4
IMMUNOGEN	Synthesized peptide derived from human CREB-2 around the non-phosphorylation site of S245.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.