

CREB3L2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11230	SYNONYMS
SIZE	100 ug	CREB3L2; BBF2H7; Cyclic AMP-responsive element-binding protein 3-like protein 2; cAMP-responsive element-binding protein 3-like protein 2; BBF2 human homolog on chromosome 7
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CREB3L2
PROTEIN NAME	Cyclic AMP-responsive element-binding protein 3-like protein 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human CREB3L2. at AA range: 240-320

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.