

EKLF/CKLF/UKLF Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX11432	SYNONYMS
SIZE	100 ug	KLF1; EKLF; Krueppel-like factor 1; Erythroid
HOST	Rabbit	krueppel-like transcription factor; EKLF; KLF5; BTEB2; CKLF; IKLF; Krueppel-like factor 5; Basic transcription element-binding protein 2; BTE-binding protein 2; Colon krueppel-like factor; GC-bo

PRODUCT INFORMATION

GENE NAME	KLF1/KLF5/KLF7
PROTEIN NAME	Krueppel-like factor 1/Krueppel-like factor 5/Krueppel-like factor 7
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human EKLF/CKLF/UKLF. at AA range: 280-360

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.