

SMF Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX14466	SYNONYMS
SIZE	100 ug	HMGXB3; KIAA0194; SMF; HMG domain-containing
HOST	Rabbit	protein 3; HMG box-containing protein 3; Protein SMF

PRODUCT INFORMATION

GENE NAME	HMGXB3
PROTEIN NAME	HMG domain-containing protein 3
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human SMF. at AA range: 100-180

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.