

Olfactory receptor 51I2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15232	SYNONYMS
SIZE	100 ug	OR51I2; Olfactory receptor 51I2; Odorant receptor
HOST	Rabbit	HOR5 ^{beta} 12; Olfactory receptor OR11-38

PRODUCT INFORMATION

GENE NAME	OR51I2
PROTEIN NAME	Olfactory receptor 51I2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Olfactory receptor 51I2. at AA range: 180-260

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.