

Olfactory receptor 8K1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15241	SYNONYMS
SIZE	100 ug	OR8K1; Olfactory receptor 8K1; Olfactory receptor
HOST	Rabbit	OR11-182

PRODUCT INFORMATION

GENE NAME	OR8K1
PROTEIN NAME	Olfactory receptor 8K1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Olfactory receptor 8K1. at AA range: 180-260

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	IF,ELISA
RECOMMENDED DILUTION	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.