

Laminin-R Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15269	SYNONYMS
SIZE	100 ug	RPSA; LAMBR; LAMR1; 40S ribosomal protein SA; 37
HOST	Rabbit	kDa laminin receptor precursor; 37LRP; 37/67 kDa laminin receptor; LRP/LR; 67 kDa laminin receptor; 67LR; Colon carcinoma laminin-binding protein; Laminin receptor 1; LamR; Laminin-binding pr

PRODUCT INFORMATION

GENE NAME	RPSA
PROTEIN NAME	40S ribosomal protein SA
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Laminin-R.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.