

## Lck Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX15279	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	LCK; Tyrosine-protein kinase Lck; Leukocyte C-terminal Src kinase; LSK; Lymphocyte cell-specific protein-tyrosine kinase; Protein YT16; Proto-oncogene
<b>HOST</b>	Rabbit	Lck; T cell-specific protein-tyrosine kinase; p56-LCK

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	LCK
<b>PROTEIN NAME</b>	Tyrosine-protein kinase Lck
<b>IMMUNOGEN</b>	Synthesized peptide derived from human Lck around the non-phosphorylation site of Y393.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse,Rat
<b>TESTED APPLICATIONS</b>	WB,IF,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**