

Galectin-7 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15290	SYNONYMS
SIZE	100 ug	LGALS7; PIG1; LGALS7B; Galectin-7; Gal-7; HKL-14;
HOST	Rabbit	PI7; p53-induced gene 1 protein

PRODUCT INFORMATION

GENE NAME	LGALS7
PROTEIN NAME	Galectin-7
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Galectin-7. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.