

Phospho-Arrestin- β -1 (S412) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15329	SYNONYMS
SIZE	100 ug	ARRB1; ARR1; Beta-arrestin-1; Arrestin beta-1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	ARRB1
PROTEIN NAME	Beta-arrestin-1
IMMUNOGEN	Synthesized peptide derived from human Arrestin- β -1 around the phosphorylation site of S412.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.