

Smad4 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15344	SYNONYMS
SIZE	100 ug	SMAD4; DPC4; MADH4; Mothers against
HOST	Rabbit	decapentaplegic homolog 4; MAD homolog 4; Mothers against DPP homolog 4; Deletion target in pancreatic carcinoma 4; SMAD family member 4; SMAD 4; Smad4; hSMAD4

PRODUCT INFORMATION

GENE NAME	SMAD4
PROTEIN NAME	Mothers against decapentaplegic homolog 4
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human Smad4. at AA range: 40-120

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.