

Arylsulfatase D Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15374	SYNONYMS
SIZE	100 ug	ARSD; Arylsulfatase D; ASD
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	ARSD
PROTEIN NAME	Arylsulfatase D
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Arylsulfatase D. at AA range: 300-380

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	IHC,IF,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.