

## Mcl-1 Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX15388	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	MCL1; BCL2L3; Induced myeloid leukemia cell
<b>HOST</b>	Rabbit	differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	MCL1
<b>PROTEIN NAME</b>	Induced myeloid leukemia cell differentiation protein Mcl-1
<b>IMMUNOGEN</b>	Synthesized peptide derived from human Mcl-1 around the non-phosphorylation site of S159.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse,Monkey
<b>TESTED APPLICATIONS</b>	WB,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**