

Phospho-ATP-citrate synthase (S455) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15524	SYNONYMS
SIZE	100 ug	ACLY; ATP-citrate synthase; ATP-citrate; pro-S-)-lyase;
HOST	Rabbit	ACL; Citrate cleavage enzyme

PRODUCT INFORMATION

GENE NAME	ACLY
PROTEIN NAME	ATP-citrate synthase
IMMUNOGEN	Synthesized peptide derived from human ATP-citrate synthase around the phosphorylation site of S455.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.