

Cleaved-Notch 4 (V1432) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15591	SYNONYMS
SIZE	100 ug	NOTCH4; INT3; Neurogenic locus notch homolog
HOST	Rabbit	protein 4; Notch 4; hNotch4

PRODUCT INFORMATION

GENE NAME	NOTCH4
PROTEIN NAME	Neurogenic locus notch homolog protein 4
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Notch 4. at AA range: 1370-1450

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.