

AMPK γ 2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15705	SYNONYMS
SIZE	100 ug	PRKAG2; 5'-AMP-activated protein kinase subunit gamma-2; AMPK gamma2; AMPK subunit gamma-2; H91620p
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PRKAG2
PROTEIN NAME	5'-AMP-activated protein kinase subunit gamma-2
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human AMPK γ 2. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.