

RNF138 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15707	SYNONYMS
SIZE	100 ug	RNF138; NARF; HSD-4; HSD4; E3 ubiquitin-protein
HOST	Rabbit	ligase RNF138; Nemo-like kinase-associated RING finger protein; NLK-associated RING finger protein; hNARF; RING finger protein 138

PRODUCT INFORMATION

GENE NAME	RNF138
PROTEIN NAME	E3 ubiquitin-protein ligase RNF138
IMMUNOGEN	Synthesized peptide derived from the Internal region of human RNF138. at AA angle: 1-80

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.