

## Phospho-NIPA (S354) Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX15709	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	ZC3HC1; NIPA; HSPC216; Nuclear-interacting partner
<b>HOST</b>	Rabbit	of ALK; Nuclear-interacting partner of anaplastic lymphoma kinase; hNIPA; Zinc finger C3HC-type protein 1

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	ZC3HC1
<b>PROTEIN NAME</b>	Nuclear-interacting partner of ALK
<b>IMMUNOGEN</b>	Synthesized peptide derived from human NIPA around the phosphorylation site of S354.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat, Monkey
<b>TESTED APPLICATIONS</b>	WB,IF,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**