

PDGF-B Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15712	SYNONYMS
SIZE	100 ug	PDGFB; PDGF2; SIS; Platelet-derived growth factor subunit B; PDGF subunit B; PDGF-2; Platelet-derived growth factor B chain; Platelet-derived growth factor beta polypeptide; Proto-oncogene c-Sis; Becaplermin
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PDGFB
PROTEIN NAME	Platelet-derived growth factor subunit B
IMMUNOGEN	Synthesized peptide derived from the Internal region of human PDGF-B. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.