

Phospho-PDGFR- α (Y754) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15715	SYNONYMS
SIZE	100 ug	PDGFRA; PDGFR2; RHEPDGFRA; Platelet-derived growth factor receptor alpha; PDGF-R-alpha; PDGFR-alpha; Alpha platelet-derived growth factor receptor; Alpha-type platelet-derived growth factor receptor; CD140 antigen-like family member A; CD14
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PDGFRA
PROTEIN NAME	Platelet-derived growth factor receptor alpha
IMMUNOGEN	Synthesized peptide derived from human PDGFR- α around the phosphorylation site of Y754.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.