

Phospho-HDAC7 (S155) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15720	SYNONYMS
SIZE	100 ug	HDAC7; HDAC7A; Histone deacetylase 7; HD7;
HOST	Rabbit	Histone deacetylase 7A; HD7a

PRODUCT INFORMATION

GENE NAME	HDAC7
PROTEIN NAME	Histone deacetylase 7
IMMUNOGEN	Synthesized peptide derived from human HDAC7 around the phosphorylation site of S155.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.