

ARMCX3 Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|--|
| CATALOG NO. | PTX15721 | SYNONYMS |
| SIZE | 100 ug | ARMCX3; ALEX3; BM-017; Armadillo repeat-containing X-linked protein 3; ARM protein lost in epithelial cancers on chromosome X 3; Protein ALEX3 |
| HOST | Rabbit | |

PRODUCT INFORMATION

| | |
|---------------------|--|
| GENE NAME | ARMCX3 |
| PROTEIN NAME | Armadillo repeat-containing X-linked protein 3 |
| IMMUNOGEN | Synthesized peptide derived from the C-terminal region of human ARMCX3. at AA range: 260-340 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human,Mouse,Rat |
| TESTED APPLICATIONS | WB,IHC,IF,ELISA |
| RECOMMENDED DILUTION | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.